

	Type	L #	Hits	Search Text	DBs	Time Stamp	Comments	Errors
1	BRS	L 1	2894	rifampin or rifaldazine or rifamycin or rifampicin	USPA T; EPO; JPO; DERW ENT	2003/06/2 9 04:46		0
2	BRS	L 6	1329 2	(zinc adj (sulphate\$1 or sulfate\$1)) or ZnSO?sub.4 or (Zn adj SO?sub.4) or ZnSO4	USPA T; EPO; JPO; DERW ENT	2003/06/2 9 04:46		0
3	BRS	L 1 1	409	(acetonitrile or propanol or isopropanol or benzene or toluene or dichloromethane or chloroform) same 6	USPA T; EPO; JPO; DERW ENT	2003/06/2 9 04:49		0
4	BRS	L 1 6	16	1 and 11	USPA T; EPO; JPO; DERW ENT	2003/06/2 9 04:48		0
5	BRS	L 2 1	69	11 and ((protein\$1 near5 (precipitat\$5 or salt\$3 or remov\$3)) or deprotein\$8)	USPA T; US-P GPUB ; EPO; JPO; DERW ENT	2003/06/2 9 05:03		0

	Type	L #	Hits	Search Text	DBs	Time Stamp	Com men ts	Er ro rs
6	BRS	L 2 7	49	21 and (hydrophobic or drug or antibiotic)	USPA T; US-P GPUB ; EPO; JPO; DERW ENT	2003/06/2 9 05:02		0
7	BRS	L 3 3	3	5135875.URPN.	USPA T	2003/06/2 9 05:02		0
8	BRS	L 3 4	5	("3949072"   "4022880"   "4160821"   "4339432"   "4734378").PN.	USPA T	2003/06/2 9 05:02		0
9	BRS	L 3 5	5	(33 or 34) and 6	USPA T; US-P GPUB ; EPO; JPO; DERW ENT	2003/06/2 9 05:02		0
10	BRS	L 4 1	2	35 and ((protein\$1 near5 (precipitat\$5 or salt\$3 or remov\$3)) or deprotein\$8)	USPA T; US-P GPUB ; EPO; JPO; DERW ENT	2003/06/2 9 05:03		0

L3 ANSWER 1 OF 1 CAPLUS COPYRIGHT 2003 ACS  
AN 1997:56454 CAPLUS  
DN 126:109002  
TI Study on HPLC assay for the plasma concentration of **rifampicin**  
and its **pharmacokinetics** of **microsphere** formulation  
AU Zhang, Wanguo; Jiang, Xuetao; Zhu, Caijuan  
CS College of Pharmacy, The Second Military Medical University, Shanghai,  
200433, Peop. Rep. China  
SO Zhongguo Kangshengsu Zazhi (1996), 21(4), 273-276  
CODEN: ZKZAEY; ISSN: 1001-8689  
PB Zhongguo Kangshengsu Zazhishe  
DT Journal  
LA Chinese  
AB Rifampicin (RFP) concn. in rabbit plasma was detd. by a HPLC assay. The  
recovery rate was 102.04%, the linear range was 0.4-12.0 .mu.g mL-1 and  
the RSD was <2%. The pharmacokinetics study of RFP in rabbits showed that  
microsphere formulation could sustain the drug release, thus injection of  
REP microspheres in rabbits gave a more stable and long-lasting plasma  
concn.

L6 ANSWER 7 OF 9 CAPLUS COPYRIGHT 2001 ACS

AN 1992:15243 CAPLUS

DN 116:15243

TI Reversed-phase liquid chromatographic method for the simultaneous determination of the antimalarial drugs sulfadoxine, pyrimethamine, mefloquine and its major carboxylic metabolite in plasma

AU Bergqvist, Yngve; Eckerbom, Solveig; Larsson, Helena; Malekzadeh, Monireh  
CS Dep. Clin. Chem., Falun Cent. Hosp., Falun, S-791 82, Swed.

SO J. Chromatogr. (1991), 571(1-2), 169-77

CODEN: JOCRAM; ISSN: 0021-9673

DT Journal

LA English

AB A high-performance liq. chromatog. method for the simultaneous detn. of sulfadoxine, pyrimethamine, mefloquine and the carboxylic metabolite of mefloquine in plasma is described. After the **proteins** have been **pptd.** with a combination of **zinc sulfate** and **acetonitrile** contg. two internal stds., pyrimethamine and mefloquine are extd. as bases and sulfadoxine and the carboxylic metabolite of mefloquine as ion-pairs with tetrabutylammonium. The drugs are sepd. by HPLC on a 3 .mu.m octadecylsilica column with UV detection at 229 nm. The method is simple and reliable and enables the simultaneous detn. of the drugs in 600-.mu.L plasma samples with a sensitivity suitable for std. drug monitoring purposes.

L6 ANSWER 9 OF 9 CAPLUS COPYRIGHT 2001 ACS

AN 1990:30169 CAPLUS

DN 112:30169

TI An improved micro-scale **protein precipitation**  
procedure for HPLC assay of therapeutic drugs in serum

AU Lam, Stanley; Malikin, Galina

CS Albert Einstein Coll. Med., Bronx, NY, 10461, USA

SO J. Liq. Chromatogr. (1989), 12(10), 1851-72

CODEN: JLCHD8; ISSN: 0148-3919

DT Journal

LA English

AB A **protein pptn.** procedure for prepg. serum-free supernatant for HPLC of therapeutic drugs is described. **Protein pptn.** is facilitated by adding small amts. of **zinc sulfate** to the serum followed by a polar org. solvent (methanol, **acetonitrile**) with subsequent centrifugation. Since the procedure involves few pipetting steps, sample loss is minimized and recovery and precision are improved. Correlation coeffs. of 1-5% are accomplished for the assays without internal stds. The **protein pptn.** procedure is applicable to the HPLC of drugs in serum at min. detection levels of 0.5 mg/mL and 0.1 mg/mL by UV and fluorescence detection, resp. The method has been applied to the detn. of several drugs in human blood serum.